


# MFAD USER MANUAL


Read this manual before operating the device



## Preparation

Wait for the chamber of the device to reach 50 °C–60 °C; the device will maintain the chamber's temperature at this range. The device will not start the disinfection cycle if the mentioned condition is not met.

- 1.) Open the door then insert the facemasks inside the chamber and hang it on the hooks located at the roof of the chamber
- 2.) Close the door and push the button to start the disinfection cycle. The UV-C will turn on and the countdown of 30 minutes will start.
- 3.) After the 30 minute disinfection cycle, the buzzer will alarm indicating that you can now retrieve the facemasks.



**Caution: Keep the device out of reach of children**

The LCD screen at the left of the door above the button will show the countdown timer